Q:\Ver10Gen\Aaadocs\Software\File_Formats\ImageIPR.doc

Image Properties (.IPR) File Format Version 3.34 July 2006				
	, ,	•		
	<u>DATA</u>	OFFSET [BYTES]	DATA <u>TYPE</u>	<u>DESCRIPTION</u>
	Version	0	short	334
	Image type	2	short	0=empty;1=electron;2=xmap;3=disk;
	5 11			4=overlay
	Label	4	char[8]	Image label
	Smin	12	short	Min collected signal
	Smax	14	short	Max collected signal
	Color	16	short	color:0=gray;1=R;2=G;3=B;4=Y5=
				M;6=C;8=overlay
	Preset mode	18	short	0=clock;1=live
	Preset time	20	long	Dwell time for x-ray (millisec)
	Data type	24	short	0=ROI; 1=Net inte;
	• •			2=K ratio; 3=Wt%; 4=Mthin
	Time constant-old	26	short	Amplifier pulse processing time [usec]
	Reserved	28	short	Not used
	ROI start channel	30	short	ROI starting channel
	ROI end channel	32	short	ROI ending channel
	User Set Minimum	34	Short	User Defined Min signal range
	User Set Maximum	36	Short	User Defined Max signal range
	iADC	38	short	Electron detector number: 1;2;3;4
	Reserved	40	short	Not used
	iBits	42	short	conversion type: 8;12 (not used)
	nReads	44	short	No. of reads per point
	nFrames	46	short	No. of frames averaged (not used)
	fDwell	48	float	Dwell time (not used)
	Accelerating Voltage	52	short	In units of 100V
	Tilt	54	short	Sample tilt [deg]
	Tkoff	56	short	Takeoff angle [deg]
	Magn	58	long	Magnification
	WD	62	short	Working distance [mm]
	mppX	64	float	Microns per pixel in X direction
	mppY	68	float	Microns per pixel in Y direction
	nTextLines	72	short	No. of comment lines
	Char text	74	char[4][32]	Comment text
	Reserved	202	float[4]	Not used
	nOverlayElements	218	short	No. of overlay elements
	0 1 1	220	1 4(1/2)	Λ 1 1 ⁻¹

Total Size: 264 bytes.

220

252

256

Overlay colors

Reserved

TimeConstant - new

short[16]

float

float[2]

Overlay colors

Not used

Amplifier time constant [usec]

Q:\Ver21\Aaadocs\Software\File_Formats\imgipr32.doc

Image Properties (.IPR) File Format Version 3.0 29-JUL-98

DATA OFFSET [BYTES] DATA TYPE **DESCRIPTION** Edit Log

Variable	Version	Purpose
User Set Minimum User Set Maximum	V334 - July 2006	To allow users to define the min/max video intensity for display. Supported in SpecUTil.
Accelerating Voltage	V333 - Sep 19, 2001	Units changed from KV to 100V
	V200 - 1998-2001	